

# High Energy Resolution Monochromatic XPS Spectra of Polyethylene

## Key Words

- Surface Analysis
- High Resolution
- High Sensitivity

Thermo Scientific ESCALAB 250 with its microfocused monochromator and coaxial charge compensation combines high energy resolution with high sensitivity. This combination is particularly important in the analysis of polymers where there is the risk of damage to the sample surface during analysis. The peak fitting routine on the *Avantage* data system is used to reveal the vibrational fine structure in the spectrum.

The high sensitivity lens allows the analysis of weaker features in the spectrum with short acquisition times. This means that the analyst can now make use of valence band spectra to distinguish very similar polymeric materials.

ESCALAB 250 utilizes a 500 mm Rowland circle monochromator to produce a small X-ray spot with minimum energy spread. When combined with the lens and coaxial charge compensation system, the very highest energy resolution and sensitivity may be obtained from polymers and other insulating materials.

The C 1s spectrum shown was obtained from clean polyethylene film, using the microfocused monochromator. The total acquisition time for the spectrum was 40 seconds. The FWHM of the total peak is 0.79 eV, but the peak shows asymmetry attributed to vibrational fine structure (Figure 1). Peak fitting has been used to resolve the separate components which each have a FWHM of 0.58 eV.

Valence band spectra can be used to 'fingerprint' otherwise similar polymers. The valence band is a relatively weak feature and traditionally acquisition times have extended over several hours. The Figure 2 spectrum shows the valence band from polyethylene, acquired in only 3 minutes 20 seconds.




Figure 1: C 1s spectrum of Polyethylene




Figure 2: Valence band spectrum of Polyethylene

## Reference

G. Beamson, D.T. Clark, J. Kendrick, D. Briggs, *J. Elec. Spec. & Rel. Phenom.* 57 (1991) 79

[www.thermo.com/surfaceanalysis](http://www.thermo.com/surfaceanalysis)

*In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.*

**Africa**  
+43 1 333 5034 127

**Australia**  
+61 2 8844 9500

**Austria**  
+43 1 333 50340

**Belgium**  
+32 2 482 30 30

**Canada**  
+1 800 530 8447

**China**  
+86 10 8419 3588

**Denmark**  
+45 70 23 62 60

**Europe-Other**  
+43 1 333 5034 127

**France**  
+33 1 60 92 48 00

**Germany**  
+49 6103 408 1014

**India**  
+91 22 6742 9434

**Italy**  
+39 02 950 591

**Japan**  
+81 45 453 9100

**Latin America**  
+1 608 276 5659

**Middle East**  
+43 1 333 5034 127

**Netherlands**  
+31 76 579 55 55

**South Africa**  
+27 11 570 1840

**Spain**  
+34 914 845 965

**Sweden/Norway/  
Finland**  
+46 8 556 468 00

**Switzerland**  
+41 61 48784 00

**UK**  
+44 1442 233555

**USA**  
+1 800 532 4752

**[www.thermo.com](http://www.thermo.com)**



*VG Systems Ltd. Trading as  
Thermo Fisher Scientific, East  
Grinstead, UK is ISO Certified.*

AN31043\_E 05/08M

*©2004, 2008 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.  
Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.*